

DNBSEQ™ SERVICE OVERVIEW

Long Non-coding RNA Sequencing



Service Description

Long non-coding RNAs (lncRNAs) are a large class of transcribed RNA molecules with a length of more than 200 nucleotides that do not encode proteins. lncRNAs are thought to encompass nearly 30,000 different transcripts in humans, hence lncRNA transcripts account for the major part of the non-coding transcriptome.

lncRNA discovery is still at an early stage and only a small proportion of lncRNAs have so far been investigated. Although we can start to classify different types of lncRNA functions, we are still far from being able to predict the function of new lncRNAs.

Sequencing Service Description

We offer expression profiling as one way to uncover the function of lncRNA. Identifying lncRNAs that are differentially expressed during development or in particular situations can shed light on their potential functions. Alternatively, looking for lncRNAs and protein-coding genes whose expression is correlated, can indicate co-regulation or related functions.

Sequencing Service Specification

DNBSEQ™ Long non-coding RNA Sequencing services are executed with the DNBSEQ™ Technology Platform, with the benefit of lower cost and faster turnaround time.



- Paired-end 100 bp
- Strand specific library with rRNA removal
- Data and bioinformatics analysis are available in standard file formats
- Advanced RNA data visualization and data mining with Dr. Tom system



- Sequencing Quality Standard
- Guaranteed ≥85% of bases with quality score of Q30



- Turnaround Time
- Typical 27 working days from sample QC acceptance to filtered data availability
 - Expedited services are available, contact our specialist for details

Project Workflow

We care for your project from the receipt of samples through to the reporting of results. Highly experienced laboratory professionals follow strict quality procedures to ensure the integrity of your results.

SAMPLE PREPARATION

Sample QC

LIBRARY CONSTRUCTION

Library QC

SEQUENCING

Sequencing QC

DATA OUTPUT

Data QC

BIOINFORMATICS ANALYSIS

Delivery QC



Quality Data



Fast TAT



Cost Effective



Data Analysis

In addition to data output, we offer a range of standard and customized bioinformatics pipelines for your Long non-coding RNA sequencing project. Reports and output data files are delivered in industry standard file formats: FASTQ, BAM and Excel.

STANDARD ANALYSIS

- Identification of mRNA and LncRNA
- Quantification and differential expression analysis
- Gene ontology analysis and pathway enrichment analysis

DR. TOM SYSTEM ANALYSIS

- Transcription factor prediction (For AnimalTFDB/PlantTFDB)
- GSEA Analysis
- Rfam、Pfam、Reactome、COG、EggNOG and InterPro annotation
- miRNA-mRNA interaction, lncRNA-mRNA interaction analysis
- Protein-Protein Interaction (PPI) analysis
- Co-expression Interaction Network Analysis
- Custom dataset table upload

Sample Requirements

We can process your total RNA, blood, cell line, FFPE, fresh frozen tissues samples from a variety of species, with the following general requirements.

Purity: OD260/280 \geq 2.0, OD260/230 \geq 2.0, without degradation and DNA contamination

Species	RNA Amount and Concentration	Quantitative Result
Human, mouse and rat	Total RNA \geq 200 ng	RIN \geq 6.0
	Concentration \geq 20 ng/ μ L	28S/18S \geq 1.0

Request for Information or Quotation

Contact a BGI Genomics representative to discuss how we can meet your specific needs or for expert advice on experiment design, from sample to bioinformatics.

info@bgi.com

www.bgi.com

For Research Use Only. Not for use in diagnostic procedures (except as specifically noted).

Copyright© BGI Genomics 2024. All trademarks are the property of BGI Genomics or their respective owners. This material contains information on products targeted to a wide range of audiences and could contain product details or information otherwise not accessible or valid in your country. Please be aware that we do not take any responsibility for accessing such information, which may not comply with any legal process, regulation, registration, or usage in the country of your origin. Unless otherwise informed, certain sequencers and sequencing reagents are not available in selected countries or regions. Please get in touch with a representative for regional availability. The company reserves the right of final interpretation.



www.bgi.com



BGI Genomics

